




NEW ENDODONTIC INSTRUMENTS

Micro Mirrors (Rhodium)




Endo Mirror Handle



#CSHS
With grooves and gauge to measure endo files.

Designed by Dr. R. Schmidt

M.T.A. Carriers



REPLACEMENT NEEDLES

#IDOVV (1.6mm) #IDOVY (.99mm)
#IDOVB (.8 mm)

REPLACEMENT RINGS

#IDOVVP (1.6mm) #IDOVRY (.99mm)
#IDOVRB (.8 mm)

REPLACEMENT BODY

#IDOVMB (ALL)

Designed by Dr. J. Dovgan

ENDODONTIC FORCEPS



Designed by Dr. M. Huey

M.T.A. Carriers (Block & Carver)



Designed by Dr. Ed Lee

2372 Stanwell Circle
P.O. Box 5988
Concord, CA 94520 USA

Phone (925) 798-2206
(800) 950-2206
Fax (925) 798-2053