



DANVILLE MATERIALS

**INNOVATIVE DENTAL
PRODUCTS**

**OTHER INNOVATIVE
PRODUCTS FROM
DANVILLE ENGINEERING**



**MICROCAB & MACROCAB
DUST CABINETS**

MICROETCHERS & ACCESSORIES

**MICROPROPHY AIR POLISHERS
& ACCESSORIES**

ABRASIVES

**MICROETCHER, MICROPROPHY &
PREPSTAR INSTALLATION KITS**

**DUST CONFINEMENT
CABINETS AND SYSTEMS**

**BONDING ENHANCEMENT
SYSTEMS**

**HIGH PRECISION
ULTRASONIC SCALERS**



DANVILLE MATERIALS

**DESENSITIZING & CARIES
DETECTION AGENTS**

**AIR INHIBITORS,
ADHESIVES & ETCHANTS**

**PROVISIONAL COMPOSITES &
CEMENTS**

FLOWABLE COMPOSITES

IMPRESSION MATERIALS

**SECTIONAL MATRIX SYSTEMS
& BITINE RINGS**

POSTS

LUTING CEMENTS

**ORTHODONTIC BONDING
MATERIALS**

AIR ABRASION

PREPSTART



**THE NON-TRAUMATIC ALTERNATIVE TO THE
HIGH-SPEED HANDPIECE.**

**Highly Rated
by
Major
Reviewers!**

The PrepStart™ air abrasion system, offers all the performance and controllability of more expensive systems for as little as one-quarter the cost. So its not surprising that PrepStart payback as a practice builder takes only a few weeks. Many dentists report a 30% increase in business with the new level of patient trust and acceptance made possible by the introduction of air abrasion. PrepStart's price is not only right, so is its size. Weighing in at only 4-1/2 pounds and taking up less than half a square foot of counter space, PrepStart could fit right on your instrument tray.

- Ultra compact.
- Crisp ON/OFF control
- 3-Stage Air Filtration
- Adjustable Abrasive Flow
- Adjustable Pressure Control
- 27u or 50u Aluminum Oxide

Portable enough to move from operatory to operatory, PrepStart connects to existing air supplies. But since it practically pays for itself in just a few weeks, you may want a PrepStart for each operatory. PrepStart can cut preps as small as half the width of your smallest bur, even incipient lesions can now be conservatively treated in seconds. You'll not only be able to treat decay at its earliest stages, you'll be able to handle procedures previously untreatable. Virtually eliminating the time and discomfort of the dental injection. PrepStart procedures take only a matter of minutes or seconds, chair turnover increases and patient anxiety is greatly reduced.



THE PREPSTART STARTER KIT
COMPLETE RESTORATIVE STARTER KIT
\$179.00 Value

**OUR CUSTOMER SERVICE SPECIALISTS
CAN HELP YOU WITH YOUR PRODUCT
AND TECHNIQUES QUESTIONS.
CALL 1-800-827-7940**



INNOVATIVE DENTAL PRODUCTS



■
GIVES YOUR EXISTING
IN-HOUSE AIR A BOOST
TO 135 PSI

PREPSTART ACCESSORY

The PowerPlus™ air booster is the easiest and fastest way of providing an upper range of air pressure for maximum performance. Many dentists find 80 to 90 psi to be more than adequate (as most compressors are now preset). However, 100 to 110 psi will cut faster. Normally, it is easy to increase the pressure range in your office by adjusting the controller box on the compressor. If, however, this adjustment is not possible, the PowerPlus or Prime Air is the perfect solution. **With the PowerPlus you don't have to touch your compressor!** The PowerPlus plugs into the existing Microetcher or similar air outlet (supplied with PrepStart). The PowerPlus doubles the air pressure up to 135 psi, starting with a 60 psi. No electricity is used, only compressed air. **The PrepStart PowerPlus combination allows a full range of PrepStart efficiency and performance.** It's operation is quiet, with a slow ticking sound when in use. The PowerPlus functions only when the PrepStart is used. Connection takes only seconds the PrepStart simply plugs into the PowerPlus and sits on top.

- Dimensions
PrepStart: 6.5" x 5.5" x 9.0"
PowerPlus: 6.7" x 3.5" x 8"
- Weight
PrepStart: 4.5 lbs.
PowerPlus: 3.5 lbs.
- PrepStart
Pressure Settings: 60-120 psi
Powder Size: 27 & 50 micron aluminum oxide
Powder Rate: 0.7 to 4.2 gm/min
- PowerPlus Power Requirements:
Min. Input Pressure: 60 psi
Max. Output Pressure: 135 psi
Max. Operating Pressure: 20 kgf
Max. Flow: 400 N l/min
Ambient & Fluid Temperature: 0-50° C (without freeze)



TO ORDER CALL
1-800-827-7940
OR CONTACT YOUR LOCAL DEALER

ORDERING INFORMATION

200830 PowerPlus air booster

200259 PrepStart air abrasion

Nozzle Options

186110	0° with .048 Tip
186112	60° with .048 Tip
186109	90° with .048 Tip
NozKit	Nozzle Kit
(Specify 3 of .048 Tip Options)	
186113	60° with .032 Tip
186114	90° with .032 Tip

Abrasives

Aluminum Oxide	
15301	50 Micron, 1 lb white
15302	50 Micron, 5 lb white
80042	27 Micron, 1 lb white

■
All Danville products carry a risk-free 30 day,
unconditional money back guarantee.*

■
**CONSULT OUR PRICE SHEET
FOR PRICING AND ADDITIONAL
PRODUCT LISTINGS.**

*Call for limitations.

Or Contact Your Dealer:

